



# Anti-Ovalbumin (native) polyclonal antibody (DPATB-H81858)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Ovalbumin Polyclonal Antibody
<b>Target</b>	Ovalbumin
<b>Immunogen</b>	Full length native protein (purified from Duck egg yolk).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Duck
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB, IEP, ID
<b>Cellular Localization</b>	Secreted.
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	0.14M Sodium chloride, 10mM PBS, pH 7.2
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

# BACKGROUND

Introduction	Storage protein of egg white. Lack protease inhibitory activity.
Keywords	Allergen Gal d 2; Allergen Gal d II; CAP 2; CAP2; Cytoplasmic antiproteinase 2; Egg albumin; Gal d II; Hen egg white; OVAL_CHICK; Ovalbumin; PI8; Plakalbumin; Protease inhibitor 8; Serpin B8; serpin peptidase inhibitor clade B member 8; SERPINB14; SERPINB