



Anti-Oxazepam polyclonal antibody (DPATB-H81859)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Oxazepam Polyclonal Antibody
Target	Oxazepam
Immunogen	Oxazepam(N-1) conjugated to BSA.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	150M Sodium chloride, 20mM PBS, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Oxazepam is a tranquilizer member of the benzodiazepine class. Oxazepam is prescribed to treat feelings of tension and anxiety. It is also used to calm patients who are suffering from the symptoms of alcohol withdrawal.

Keywords

Alepam; Medopam; Murelax; Noripam; Opamox; Ox-Pam; Purata; Serax; Serepax, Vaben; Sobril; Oxascand; Zaxpam ; Sheep Anti-Oxazepam Polyclonal Antibody; Anti-Oxazepam Polyclonal Antibody; Oxazepam Polyclonal Antibody Sheep Anti-Oxazepam PAb; Anti-Oxazepam PA