



Anti-P. aeruginosa type Ds Group Ds Polyclonal antibody (DPATB-H81862)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Pseudomonas aeruginosa Group Ds Polyclonal Antibody
Target	P. aeruginosa type Ds Group Ds
Immunogen	Tissue / cell preparation of Pseudomonas aeruginosa Group Ds (heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	P. aeruginosa
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	0.1% Sodium Azide
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction Pseudomonas aeruginosa is Gram-negative, aerobic, rod-shaped bacteria with unipolar motility. An opportunistic human pathogen, P. aeruginosa is also an opportunistic pathogen of plants. P. aeruginosa bacteria are clinically important because they are resi

Keywords Pseudomonas aeruginosa Group DsP aeruginosa; P. aeruginosa; Pseudomonas aeruginosa; Pseudomonadaceae; Pseudomonas; pseudomonas aeruginosa; Gamma Proteobacteri
