



Anti-*P. aeruginosa* Polyclonal antibody (DPATB-H81863)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Guinea pig Anti-Pseudomonas Polyclonal Antibody
Target	<i>P. aeruginosa</i>
Immunogen	Outer membrane protein extract of <i>P. aeruginosa</i> .
Isotype	IgG
Source/Host	Guinea pig
Species Reactivity	<i>P. aeruginosa</i>
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, Agg
Format	Liquid
Size	100 µg
Buffer	Not applicable.
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Pseudomonas aeruginosa is Gram-negative, aerobic, rod-shaped bacteria with unipolar motility. An opportunistic human pathogen, *P. aeruginosa* is also an opportunistic pathogen of plants. *P. aeruginosa* bacteria are clinically important because they are resi

Keywords

P aeruginosa; *P. aeruginosa*; *Pseudomonas aeruginosa*; *Pseudomonadaceae*; *Pseudomonas*; *pseudomonas aeruginosa*; *Gamma Proteobacteria*; *Bacteria*