



Anti-*P. aeruginosa* Polyclonal antibody (DPATB-H81864)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Pseudomonas Polyclonal Antibody
Target	<i>P. aeruginosa</i>
Immunogen	Tissue / cell preparation-whole cells of <i>Pseudomonas aeruginosa</i> .
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>P. aeruginosa</i>
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IHC-Fr, WB, ELISA
Format	Liquid
Size	1 ml
Buffer	0.01M PBS
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Pseudomonas aeruginosa is Gram-negative, aerobic, rod-shaped bacteria with unipolar motility. An opportunistic human pathogen, P. aeruginosa is also an opportunistic pathogen of plants. P. aeruginosa bacteria are clinically important because they are resi

Keywords P aeruginosa; P. aeruginosa; Pseudomonas aeruginosa; Pseudomonadaceae; Pseudomonas; pseudomonas aeruginosa; Gamma Proteobacteria; Bacteria
