



Anti-*P. aeruginosa* type Cs Group Cs Polyclonal antibody (DPATB-H81866)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Pseudomonas aeruginosa Group Cs Polyclonal Antibody
Target	<i>P. aeruginosa</i> type Cs Group Cs
Immunogen	Tissue / cell preparation of <i>Pseudomonas aeruginosa</i> Group Cs (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>P. aeruginosa</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Pseudomonas aeruginosa is Gram-negative, aerobic, rod-shaped bacteria with unipolar motility. An opportunistic human pathogen, *P. aeruginosa* is also an opportunistic pathogen of plants. *P. aeruginosa* bacteria are clinically important because they are resi

Keywords

Pseudomonas aeruginosa Group CsP *aeruginosa*; *P. aeruginosa*; *Pseudomonas aeruginosa*; *Pseudomonadaceae*; *Pseudomonas*; *pseudomonas aeruginosa*; Gamma Proteobacteri