



Anti-*P. marginalis* Polyclonal antibody (DPATB-H81870)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat <i>Pseudomonas marginalis</i> (conjugated) Polyclonal Antibody
Target	<i>P. marginalis</i>
Immunogen	Tissue / cell preparation- <i>Pseudomonas marginalis</i> total proteins
Isotype	IgG
Source/Host	Rat
Species Reactivity	<i>P. marginalis</i>
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Pseudomonas marginalis is a phytopathogenic, gram-negative bacteria, causing a disease commonly known as marginal leaf spot and also soft rot of vegetables. The presence of *P. marginalis* is an important cause of vegetable spoilage before and after harvest

Keywords

Bacteria; Proteobacteria; Gamma Proteobacteria; Pseudomonadales; *P. marginalis*; Pseudomonadaceae; *Pseudomonas*; *Pseudomonas marginalis*