



# Anti-*P. putida* Polyclonal antibody (DPATB-H81871)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit <i>Pseudomonas putida</i> (conjugated) Polyclonal Antibody
<b>Target</b>	<i>P. putida</i>
<b>Immunogen</b>	Synthetic peptide to <i>Pseudomonas putida</i> total proteins.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	<i>P. putida</i>
<b>Purification</b>	Preadsorbed antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	50% Glycerol, ddH2O
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Pseudomonas putida is a soil bacterium; is Gram-negative with rod-shaped cells and multitrichous flagella. Pseudomonas putida has the potential to help clean up organic pollutants it is involved in breaking down aromatic hydrocarbons, like TNT. Aromatic h
<b>Keywords</b>	Bacteria; Proteobacteria; Gamma Proteobacteria; Pseudomonadales; P. putida; Pseudomonadaceae; Pseudomonas; Pseudomonas putida