



Anti-*P. stutzeri* Polyclonal antibody (DPATB-H81873)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit <i>Pseudomonas stutzeri</i> (conjugated) Polyclonal Antibody
Target	<i>P. stutzeri</i>
Immunogen	Tissue / cell preparation; <i>Pseudomonas stutzeri</i> total proteins
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>P. stutzeri</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH ₂ O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Pseudomonas stutzeri is a Gram-negative, rod-shaped, motile, single polar-flagellated, soil bacterium first isolated from human spinal fluid. It is a denitrifying bacterium, and strain KC of P. stutzeri may be used for bioremediation as it is able to degr

Keywords Bacteria; Proteobacteria; Gamma Proteobacteria; Pseudomonadales; P. Stutzeri; Pseudomonadaceae; Pseudomonas; Pseudomonas stutzeri; S. aureus
