



Anti-Papillomavirus type 2 Envelope Antigen 7 Polyclonal antibody (DPATB-H81877)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Canine papillomavirus3E7 Polyclonal Antibody
Target	Papillomavirus type 2 Envelope Antigen 7
Immunogen	Synthetic peptide corresponding to C-terminal residues of Canine papillomavirus3E7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Papillomavirus
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Papillomaviridae is an ancient taxonomic family of non-enveloped DNA viruses, collectively known as papillomaviruses. Several hundred species of papillomaviruses, traditionally referred to as "types", have been identified infecting all carefully inspected

Keywords

Papillomaviridae; Canine papillomavirus3E7; E7
