



# Anti-Parathion polyclonal antibody (DPATB-H81881)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |  |
|---------------------------|--|
| <b>Product Overview</b>   | Rabbit Anti-Parathion Polyclonal Antibody  |
| <b>Target</b>             | Parathion  |
| <b>Immunogen</b>          | Chemical / Small Molecule: Parathion conjugated to KLH.  |
| <b>Isotype</b>            | IgG  |
| <b>Source/Host</b>        | Rabbit   |
| <b>Species Reactivity</b> | N/A  |
| <b>Purification</b>       | Whole antiserum  |
| <b>Conjugate</b>          | Unconjugated   |
| <b>Applications</b>       | ELISA  |
| <b>Format</b>             | Liquid   |
| <b>Size</b>               | 200 µl   |
| <b>Buffer</b>             | 25% Glycerol, Whole serum.   |
| <b>Preservative</b>       | None   |
| <b>Storage</b>            | Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles. |

## BACKGROUND

|                     |  |
|---------------------|--|
| <b>Introduction</b> | A highly poisonous yellowish liquid organic phosphate insecticide (C <sub>8</sub> H <sub>10</sub> NO <sub>5</sub> PS) that is an irreversible inhibitor of cholinesterase. |
| <b>Keywords</b>     | Parathion  |