



Anti-Phage Phage Tail Tape Measure Protein Polyclonal antibody (DPATB-H81887)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Phage Tail Tape Measure Protein Polyclonal Antibody
Target	Phage Phage Tail Tape Measure Protein
Immunogen	Synthetic peptide corresponding to N terminal residues of Staphylococcus phage phi13 Phage Tail Tape Measure Protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Phage
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Phage Tail Tape Measure Protein contains domains of phage tail tape measure protein of TP901 family, Domain of unknown function DUF20, Phage-related protein, as well as Lytic Transglycosylase (LT) and Goose Egg White Lysozyme (GEWL) domain. Members includ

Keywords

Phage Tail Tape Measure Protein
