



Anti-Quinolone polyclonal antibody (DPATB-H81936)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Quinolone Polyclonal Antibody
Target	Quinolone
Immunogen	Norfloxacin(N)conjugated to BSA
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Quinolones and fluoroquinolones form a group of broad-spectrum antibiotics. They are related to nalidixic acid. Fluoroquinolone antibiotics are highly potent and at one time they were considered relatively safe. However, they can have potentially troubles
Keywords	Flumequine; 7-fluoro-12-methyl-4-oxo-1-azatricyclo[7.3.1.0 ^{5,13}]trideca-2,5,7,9-tetraene-3-carboxylic acid; Apurone; R 802; Imequyl; Quinolone