



# Anti-serum polyclonal antibody (DPATB-H81939)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat Anti-Rabbit serum Polyclonal Antibody
<b>Target</b>	Rabbit serum
<b>Immunogen</b>	Normal rabbit serum
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Rabbit
<b>Purification</b>	IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IEP
<b>Format</b>	Lyophilised: To one vial of lyophilized powder add 3ml of deionized water.
<b>Size</b>	2 ml
<b>Buffer</b>	0.1M PBS, pH 7.2
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C. Aliquot and store at -20°C or -80°C when reconstituted. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Serum is the clear, straw-colored, liquid portion of blood plasma that does not contain fibrinogen or blood cells and remains fluid after clotting.

---

**Keywords**

Rabbit serum

---