



Anti-Rattin (aa 20-35) polyclonal antibody (DPATB-H81948)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Rattin Polyclonal Antibody
Target	Rattin
Immunogen	Synthetic peptide corresponding to internal amino acids 20-35 (VKRLESPNKTRRPYGA) of Rat Rattin peptide.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted
Positive Control	Rattin peptide or recombinant Rattin peptide produced in transfected 293T cells
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	0.02% Sodium Azide

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Rattin is a novel rat gene which encodes a peptide homologous to Humanin, a secreted peptide that specifically protects against neuronal cell death induced by beta-amyloid peptide (Ab) or by mutations causing early-onset familial Alzheimer's disease. Ratt

Keywords

Rattin
