



# Rabbit Anti-Reovirus p17 Ppolyclonal antibody (DPATB-H81951)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-Reovirus p17 protein Polyclonal Antibody
Target	Reovirus p17 protein
Immunogen	Synthetic peptide of Reovirus p17 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Virus
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	100 µg
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** The biological role of Reovirus p17 protein is fully unknown so far. It is possible that Reovirus

p17 protein possesses a growth retardation function.

---

**Keywords**

Reovirus p17 protein; Reovirus; p17 protein

---