



# Magic™ Anti-Reticulo endotheliosis polyclonal antibody (DPATB-H81952)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Chicken Anti-Reticulo endotheliosis Polyclonal Antibody
<b>Target</b>	Reticulo endotheliosis
<b>Immunogen</b>	Reticulo endotheliosis.
<b>Isotype</b>	IgY
<b>Source/Host</b>	Chicken
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, Neut
<b>Format</b>	Liquid
<b>Size</b>	200 µl
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

The reticulo endotheliosis virus is a type of gammaretrovirus that may cause leukemia.

---

**Keywords**

Reticulo endotheliosis

---