



Anti-Spectinomycin polyclonal antibody (DPATB-H82154)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Spectinomycin Polyclonal Antibody
Target	Spectinomycin
Immunogen	Spectinomycin conjugated to BSA
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	0.33% Sodium phosphate, 0.88% Sodium chloride
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Spectinomycin

Keywords

Spectinomycin