



Anti-Spermidine polyclonal antibody (DPATB-H82155)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Spermidine Polyclonal Antibody
Target	Spermidine
Immunogen	Spermidine conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	200 µl
Buffer	Whole Serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Spermidine is a polyamine compound was originally isolated from semen. It is found ribosomes and living tissues and has various metabolic functions.

Keywords

Spermidine
