



Anti-DNase polyclonal antibody (DPATB-H82159)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Staphylococcus DNase Polyclonal Antibody
Target	Staphylococcus DNase
Immunogen	Highly purified (> 95%) Staphylococcus aureus DNase
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Staphylococcus aureus
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	250 µl
Buffer	1% BSA, 0.01M PBS, pH 7.2
Preservative	None
Storage	Aliquot and store at -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	One of the important virulence factors in <i>S. aureus</i> is DNase. This can be used to distinguish <i>Staphylococcus aureus</i> from the coagulase-negative staphylococci by the so called DNase test.
Keywords	<i>Staphylococcus</i> DNase; <i>Staphylococcus aureus</i> ; <i>S. aureus</i> Bacteria; Eubacteria; Firmicutes; Bacilli; Bacillales; Stap