



Anti-Substance P (full length) polyclonal antibody (DPATB-H82166)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Substance P Polyclonal Antibody
Target	Substance P
Immunogen	Recombinant full length protein conjugated to BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Guinea pig
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, Dot
Cellular Localization	Secreted.
Format	Liquid
Size	50 µl
Buffer	Liquid antiserum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Tachykinins are active peptides which excite neurons, evoke behavioral responses, are potent vasodilators and secretagogues, and contract (directly or indirectly) many smooth muscles.
Keywords	C-terminal-flanking peptide,Hs.2563,Neurokinin 1,Neurokinin 2,Neurokinin A,Neurokinin alpha,Neuromedin L,Neuropeptide gamma,Neuropeptide K,NK2,NKA,NKNA,NPK,PPT,Protachykinin 1 precursor,Substance K,SubstanceP,TAC1,TAC2,Tachykinin 1,Tachykinin 2,Tachykinin

GENE INFORMATION

Entrez Gene ID	6863
UniProt ID	P20366