



Anti-Succinylarginine polyclonal antibody (DPATB-H82167)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Succinylarginine Polyclonal Antibody
Target	Succinylarginine
Immunogen	Succinylarginine pentane spacer-e-N-linked to lysine on proprietary carrier peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ICC
Format	Liquid
Size	500 µl
Buffer	98.95% Phosphate Buffer, 1% Goat Serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction	Succinylarginine
Keywords	N2-succinylarginine; N2-succinyl-L-arginine; N~2~-SUCCINYLRARGININE; CPD-421; DB02501; (2S)-5-(diaminomethylideneamino)-2-[(4-hydroxy-4-oxobutanoyl)amino]pentanoic acid