



Anti-Sulphamerazine polyclonal antibody (DPATB-H82172)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Sulphamerazine Polyclonal Antibody
Target	Sulphamerazine
Immunogen	Sulphamerazine-BTG.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction Sulfamerazine is a sulfonamide antibacterial.

Keywords Sulphamerazine; sulfonamide antibacterial.
