



# Anti-SUMO 3 (full length) polyclonal antibody (DPATB-H82178)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Sumo 3 Polyclonal Antibody
<b>Target</b>	Sumo 3
<b>Immunogen</b>	Recombinant full length protein (Arabidopsis thaliana).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Cytoplasm.
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** Ubiquitin-like protein which can be covalently attached to target lysines either as a monomer or as a lysine-linked polymer. Does not seem to be involved in protein degradation and may function as an antagonist of ubiquitin in the degradation process. Pla

---

**Keywords** Small ubiquitin related modifier 3, Small ubiquitin-related modifier 3, SMT3 homolog 1, SMT3 suppressor of mif two 3 homolog 1, SMT3 suppressor of mif two 3 homolog 3, SMT3A, Smt3B, SMT3H1, SUMO-2, SUMO-3, SUMO3, SUMO3\_HUMAN, Ubiquitin like protein SMT3A, Ubiquitin-li

---

## GENE INFORMATION

**Entrez Gene ID** [6612](#)

---

**UniProt ID** [P55854](#)

---