



# Anti-Syntaxin 16-Drosophila Golgi/Trans-Golgi Marker (aa 200-300) polyclonal antibody (DPATB-H82179)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                                                                       |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>      | Rabbit Anti-Syntaxin 16-Drosophila Golgi/Trans-Golgi Marker Polyclonal Antibody                                                                                                                       |
| <b>Target</b>                | Syntaxin 16-Drosophila Golgi/Trans-Golgi Marker                                                                                                                                                       |
| <b>Immunogen</b>             | Synthetic peptide conjugated to KLH derived from within residues 200-300 of Fruit fly (Drosophila melanogaster) Syntaxin 16. Read Abcams proprietary immunogen policy (Peptide available as ab32339.) |
| <b>Isotype</b>               | IgG                                                                                                                                                                                                   |
| <b>Source/Host</b>           | Rabbit                                                                                                                                                                                                |
| <b>Species Reactivity</b>    | Drosophila                                                                                                                                                                                            |
| <b>Purification</b>          | Immunogen affinity purified                                                                                                                                                                           |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                                                                          |
| <b>Applications</b>          | IHC-P, WB                                                                                                                                                                                             |
| <b>Cellular Localization</b> | Cytoplasm and Golgi apparatus membrane.                                                                                                                                                               |
| <b>Positive Control</b>      | Schneider L2 (Fruit fly (Drosophila melanogaster)) Whole cell lysate.                                                                                                                                 |
| <b>Format</b>                | Liquid                                                                                                                                                                                                |
| <b>Size</b>                  | 100 µg                                                                                                                                                                                                |
| <b>Buffer</b>                | 1% BSA, PBS. pH 7.4                                                                                                                                                                                   |

|                     |                                                                                                                 |
|---------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>Preservative</b> | 0.02% Sodium Azide                                                                                              |
| <b>Storage</b>      | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |

## BACKGROUND

|                     |                                                                                           |
|---------------------|-------------------------------------------------------------------------------------------|
| <b>Introduction</b> | SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network. |
|---------------------|-------------------------------------------------------------------------------------------|

|                 |                                                      |
|-----------------|------------------------------------------------------|
| <b>Keywords</b> | hsyn16,MGC90328,Stx16,STX16_HUMAN,Syn16,Syntaxin-16, |
|-----------------|------------------------------------------------------|