



Anti-TBP like protein TLP (aa 198-215) polyclonal antibody (DPATB-H82181)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-TBP like protein TLP Polyclonal Antibody
Target	TBP like protein TLP
Immunogen	Synthetic peptide: MTRHQMAQHHAQQPHPQIY, corresponding to amino acids 198-215 of Caenorhabditis elegans TBP like protein TLP MTRHQMAQHHAQQPHPQIY Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Caenorhabditis elegans
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-Fr
Cellular Localization	Cytoplasm. Nucleus.
Format	Liquid
Size	200 µl
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze

/ thaw cycles.

BACKGROUND

Introduction

Does not bind the TATA box. Has DNA-binding ability.

Keywords

21 kDa TBP-like protein, 21 kDa TBP like protein, Second TBP of unique DNA, Second TBP of unique DNA protein, STUD, TATA box binding protein like protein 1, TATA box binding protein-related factor 2, TATA box-binding protein-like protein 1, TATA box-binding protei