



Anti-TFAM (native) polyclonal antibody (DPATB-H82187)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-TFAM Polyclonal Antibody
Target	TFAM
Immunogen	Native recombinant Mouse mtTFA protein with a C-terminal7-ljtttag.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Positive Control	Rat kidney lysate.
Format	Liquid
Size	100 µl
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at-20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

This gene encodes a mitochondrial transcription factor that is a key activator of mitochondrial transcription as well as a participant in mitochondrial genome replication. Studies in mice have demonstrated that this gene product is required to regulate th

Keywords

TFAM; transcription factor A, mitochondrial; TCF6, TCF6L2; Mitochondrial transcription factor 1; mitochondrial transcription factor A; MtTF1; mtTFA; TCF 6; TCF-6; TCF6; TCF6L1; TCF6L2; TCF6L3; TFAM; TFAM_HUMAN; Transcription factor 6; Transcription factor

GENE INFORMATION

Entrez Gene ID

[7019](#)

UniProt ID

[A8MRB2](#)
