



Anti-TFIIB (aa 195-294) polyclonal antibody (DPATB-H82189)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-TFIIB Polyclonal Antibody
Target	TFIIB
Immunogen	Recombinant fusion protein: TTQKKTTIKN VIKTKNLKRE DMETENEQNE AQPNEKNEA QPNNEKNKAP PNDEKNEAPS NEENVVQDLN KDELSDNLEN NPVKNFDTDF ENKNDGEEKY, corresponding to amino acids 195-294 of Plasmodium falciparum TFIIB TTQKKTTIKNVIKTKNLKREDMETENEQNEAQPNEKNEA QPNNEKN
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Plasmodium falciparum
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Size	100 µl
Buffer	50% Glycerol, Whole serum
Preservative	None

Storage

Shipped at 4°C. Upon delivery aliquot and store at-20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.
Keywords	General transcription factor IIB,General transcription factor TFIIIB,GTF 2B,GTF2B,OTTHUMP00000058838,RNA polymerase II transcription factor IIB,S300 II,S300-II,S300II,TF 2B,TF IIB,TF2B,TF2B_HUMAN,TFIIIB,Transcription initiation factor IIB,