



Anti-thylakoid membrane proteins polyclonal antibody (DPATB-H82193)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-thylakoid membrane proteins Polyclonal Antibody
Target	thylakoid
Immunogen	Extracted thylakoid membrane proteins from spinach chloroplasts
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Spinach
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	thylakoid membrane
Format	Liquid
Size	100 µl
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

A thylakoid is a membrane-bound compartment inside chloroplasts and cyanobacteria. They are the site of the light-dependent reactions of photosynthesis. Thylakoids consist of a thylakoid membrane surrounding a thylakoid lumen. Chloroplast thylakoids frequ

Keywords

thylakoid; thylakoid membrane proteins
