



Anti-thylakoid lumen proteins polyclonal antibody (DPATB-H82194)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-thylakoid lumen proteins Polyclonal Antibody
Target	thylakoid
Immunogen	Extracted thylakoid luminal proteins from spinach chloroplasts.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Spinach
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	thylakoid lumen
Positive Control	Spinach luminal proteins.
Format	Liquid
Size	100 µl
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

A thylakoid is a membrane-bound compartment inside chloroplasts and cyanobacteria. They are the site of the light-dependent reactions of photosynthesis. Thylakoids consist of a thylakoid membrane surrounding a thylakoid lumen. Chloroplast thylakoids frequ

Keywords

thylakoid; thylakoid lumen proteins
