



Anti-TNF Receptor I polyclonal antibody (DPATB-H82198)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-TNF Receptor I Polyclonal Antibody
Target	TNF Receptor I
Immunogen	Recombinant human sTNFRp55 expressed in Chinese Hamster Cells.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	IgG fraction
Conjugate	Unconjugated
Applications	Dot, ELISA, FC, IP, RIA, WB
Cellular Localization	Cell membrane. Secreted.
Format	Liquid
Size	500 µg
Buffer	0.15M Sodium Chloride, 0.02M Potassium Phosphate. pH 7.2
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Receptor for TNFSF2/TNF-alpha and homotrimeric TNFSF1/lymphotoxin-alpha. The adapter molecule FADD recruits caspase-8 to the activated receptor. The resulting death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiate
Keywords	CD120a,CD120a antigen,FPF,MGC19588,p55,p55 R,p55-R,p55R,p60,TBP1,TBPI,TNF R,TNF R I,TNF R1,TNF R55,TNF-R1,TNF-RI,TNFAR,TNFR 1,TNFR I,TNFR-I,TNFR1,TNFR55,TNFR60,TNFR1,TNFRSF1a,TNFRSF1a,TNR1A_HUMAN,Tumor necrosis factor alpha receptor,Tumor necrosis factor

GENE INFORMATION

Entrez Gene ID	7132
UniProt ID	P19438