



Anti-Triamcinolone polyclonal antibody (DPATB-H82205)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Triamcinolone Polyclonal Antibody
Target	Triamcinolone
Immunogen	Triamcinolone-BTG.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	99% Phosphate Buffer, 0.88% Sodium chloride
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Triamcinolone is a long-acting synthetic corticosteroid given orally, by injection, by inhalation, or as a topical ointment or cream.
Keywords	Triamcinolone; (11 β ,16 α)-9-fluoro-11,16,17,21-tetrahydroxypregna-1,4-diene-3,20-dione