



Anti-Vipera russelli venom polyclonal antibody (DPATB-H82276)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Vipera russelli venom Polyclonal Antibody
Target	Reacts with Vipera russelli.
Immunogen	Crude venom
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Vipera russelli is a venomous viper species, which is found in Asia throughout the Indian subcontinent, much of Southeast Asia, southern China and Taiwan. This species the most dangerous viperid snake and is a major cause of snakebite injury and mortality

Keywords

Vipera russelli; *Vipera russelli* venom
