



# Anti-M. pneumoniae Polyclonal antibody (DPATB-H82504)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Mycoplasma pneumoniae Polyclonal Antibody
<b>Target</b>	M. pneumoniae
<b>Immunogen</b>	Purified Mycoplasma pneumoniae
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	M. pneumoniae
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IP, ICC, IF
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Mycoplasma pneumoniae is a very small bacterium, a member of the class Mollicutes, meaning soft skin. Along with the other members of this class (Acholeplasma, Anaeroplasm, Asteroleplasma, Spiroplasma, and Ureaplasma) Mycoplasma are characterized by thei
---------------------	--

---

<b>Keywords</b>	bacterial pneumonia; M. pneumoniae;
-----------------	-------------------------------------

---