



Anti-Shigella Polyclonal antibody (DPATB-H82561)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Shigella Polyclonal Antibody
Target	Shigella
Immunogen	Mixture of <i>S. boydii</i> , <i>S. flexneri</i> and <i>S. dysenteriae</i> .
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Shigella
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-Fr, ICC/IF
Format	Liquid
Size	500 µl
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Shigella is a genus of Gram-negative, non-motile, non-spore forming rod-shaped bacteria closely related to <i>Escherichia coli</i> and <i>Salmonella</i> . The causative agent of human shigellosis,
--------------	--

Shigella also cause disease in other primates, but not in other mammals.

Keywords

S. boydii, S. dysenteriae, S. flexneri, S. sonnei, Shigella Boydii, Shigella Dysenteriae, Shigella Flexneri, Shigella Sonnei,
