



Anti-EBV EBNA1 Polyclonal antibody (DPATB-H82740)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-EBV EBNA-1 Polyclonal Antibody
Target	EBV EBNA1
Immunogen	Recombinant p18 of EBV EBNA-1.
Isotype	IgG
Source/Host	Goat
Species Reactivity	EBV
Conjugate	Unconjugated
Applications	WB, ICC, IF, ELISA
Format	Liquid
Size	1 ml
Buffer	pH: 7.20Preservative: 0.1% Sodium azideConstituent: 99% PBS
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction The Epstein–Barrvirus (EBV), also called human herpesvirus5(HHV-4), is a virus of the herpes

family, which includes herpes simplex virus 1 and 2, and is one of the most common viruses in humans. It is best known as the cause of infectious mononucleosis. It i

Keywords	HHV 4; HHV4; human herpesvirus 4;
-----------------	-----------------------------------
