



Anti-HSV type 2 Polyclonal antibody (DPATB-H82741)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-HSV2 Polyclonal Antibody
Target	HSV type 2
Immunogen	HSV2, strain G (Human).
Isotype	IgG
Source/Host	Sheep
Species Reactivity	HSV
Conjugate	Unconjugated
Applications	WB, ELISA, ICC, IF
Format	Liquid
Size	1 ml
Buffer	pH: 7.20Preservative: 0.1% Sodium azideConstituent: 99% PBS
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Herpes simplex type 3 (HSV2) belongs to a family that includes HSV1, Epstein-Barr virus (EBV)

and Varicella zoster (chicken pox) virus. HSV1 and HSV2 are extremely difficult to distinguish from each other. These viruses have a DNA genome, an icosahedral pro

Keywords

Herpes Simplex Virus Type 2;
