



Anti-HSV type 1 Polyclonal antibody (DPATB-H82742)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-HSV1 Polyclonal Antibody
Target	HSV type 1
Immunogen	HSV type 1, strain F.
Isotype	IgG
Source/Host	Goat
Species Reactivity	HSV
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ICC, IF
Format	Liquid
Size	1 ml
Buffer	pH: 7.20Preservative: 0.1% Sodium azideConstituent: 99% PBS
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Herpes simplex type 1 (HSV-1) belongs to a family that includes HSV-2, Epstein-Barr virus (EBV) and Varicella zoster (chicken pox) virus amongst others. HSV-1 and HSV-3 are extremely difficult to distinguish from each other. Members of this family have a?

Keywords

Herpes simplex virus 1; Herpes simplex virus; Herpes simplex virus type 1; HSV 1; Human herpesvirus 1; Human herpesvirus type 1;
