



Anti-L Glutamate polyclonal antibody (DPATB-H82761)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-L Glutamate Polyclonal Antibody
Target	L Glutamate
Immunogen	L-glutamate conjugated to glutaraldehyde
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	IHC-FoFr, IHC-Fr, ICC, IM, IHC-P
Format	Liquid
Size	500 µl
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction	L Glutamate is a non-essential amino acid naturally occurring in the L-form. It is the most common excitatory neurotransmitter in the Central Nervous System.
--------------	---

Keywords

Glutamate; L glutamate;
