



# Anti-PTENP1 (full length) polyclonal antibody (DPATB-H82763)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Mouse Anti-PTENP1 Polyclonal Antibody
<b>Target</b>	PTENP1
<b>Immunogen</b>	Full length PTENP1 protein (Human)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Cytoplasmic
<b>Positive Control</b>	PTENP1 transfected 293T cell lysate
<b>Format</b>	Liquid
<b>Size</b>	50 µg
<b>Buffer</b>	PBS, pH 7.2
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.

---

## BACKGROUND

### Introduction

PTENP1 is also called phosphatase and tensin homolog pseudogene 1. Magnesium is its cofactor.

---

### Keywords

Phosphatase and tensin homolog pseudogene 1; psiPTEN; PTEN rs; PTEN2; PTH2; TEP1;

---