



Anti-vjbR (full length) polyclonal antibody (DPATB-H82765)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-vjbR Polyclonal Antibody
Target	vjbR
Immunogen	Full length Brucella suis vjbR protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Brucella
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

vjbR is a quorum-sensing regulator that activates the expression of the virB secretion system and the flagellar apparatus of *Brucella melitensis*. The vjbR (HTH-type quorum sensing-dependent transcriptional regulator vjbR) functions as a transcriptional re

Keywords

HTH type quorum sensing dependent transcriptional regulator vjbR,