



# Anti-Tail Sheath (N-terminal) polyclonal antibody (DPATB-H82766)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Tail Sheath Polyclonal Antibody
<b>Target</b>	Tail Sheath
<b>Immunogen</b>	Synthetic peptide corresponding to the very N-terminal residues of Pseudomonas aeruginosa PSPA7_0697
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Pseudomonas
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Format</b>	Liquid
<b>Size</b>	50 µg
<b>Buffer</b>	50% Glycerol, PBS
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

PSPA7\_0697 is a phage tail sheath protein.

---

**Keywords**

PSPA7\_5059,Tail sheath,Tail sheath protein,

---