



Anti-Tail fiber protein (N-terminal) polyclonal antibody (DPATB-H82767)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Tail fiber protein Polyclonal Antibody
Target	Tail fiber protein
Immunogen	Synthetic peptide corresponding to residues near the N terminus of <i>Pseudomonas aeruginosa</i> Tail fiber protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Pseudomonas</i>
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Tail fiber protein is a phage-related tail fibre protein.

Keywords

PSPA7_0787,
