



Anti-Tail Sheath (C-terminal) polyclonal antibody (DPATB-H82769)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Tail Sheath Polyclonal Antibody
Target	Tail Sheath
Immunogen	Synthetic peptide corresponding to C terminal residues of bacteria <i>Pseudomonas aeruginosa</i> PSPA7_0697.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Pseudomonas</i>
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction PSPA7_0698 is a phage tail sheath protein.

Keywords PSPA7_5059,Tail sheath,Tail sheath protein,
