



Anti-C. difficile Cell Surface Protein (C-terminal) polyclonal antibody (DPATB-H82779)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Clostridium difficile Cell Surface Protein Polyclonal Antibody
Target	Clostridium difficile Cell Surface Protein
Immunogen	Synthetic peptide corresponding to C-terminal residues of Clostridium difficile Cell Surface Protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Clostridium difficile
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Clostridium difficile Cell Surface Protein contains a CW_binding_2 domain (Putative cell wall binding repeat 2) and a MurNAc-LAA domain: N-acetylmuramoyl-L-alanine amidase or MurNAc-LAA (also known as peptidoglycan aminohydrolase, NAMLA amidase, NAMLAA, A
Keywords	Cell surface protein (N acetylmuramoyl L alanine amidase);