



Anti-FIR1 polyclonal antibody (DPATB-H82784)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-FIR1 Polyclonal Antibody
Target	FIR1
Immunogen	Synthetic peptide corresponding to middle amino acid residues of <i>Pseudomonas aeruginosa</i> FIR1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Pseudomonas</i>
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction FIR1 is a phage tail sheath protein.

Keywords

Tail sheath; Tail sheath protein;
