



Anti-Helicobacter pylori Polysaccharide Deacetylase (N-terminal) polyclonal antibody (DPATB-H82789)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Helicobacter pylori Polysaccharide Deacetylase Polyclonal Antibody
Target	Helicobacter pylori Polysaccharide Deacetylase
Immunogen	Synthetic peptide corresponding to N terminal residues of Helicobacter pylori P12 Polysaccharide Deacetylase (UniProt ID: B6JKN9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Helicobacter pylori
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction	Polysaccharide Deacetylase contains CDA1 which is predicted xylanase/chitin deacetylase [Carbohydrate transport and metabolism] and domain of Polysacc deac 1; Polysaccharide deacetylase. Polysaccharide Deacetylase functions in biosynthesis of surface poly
Keywords	HPP12 0309; HPP12 0309 polysaccharide deacetylase; Polysaccharide deacetylase;