



Anti-COPV E6 (N-terminal) polyclonal antibody (DPATB-H82796)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-COPVF7 Polyclonal Antibody
Target	COPV E6
Immunogen	Synthetic peptide corresponding to the N terminal residues of COPV E6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Host Nuclear Matrix
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction

COPVF7(Canine oral papillomavirus Protein E6) is a transcriptional transactivator which binds double stranded DNA (by similarity). COPVF7is localized in host nucleus matrix. It belongs to the papillomaviridaeF7protein family.

Keywords

E6; Early protein;
