



# Anti-ULP-1 polyclonal antibody (DPATB-H82808)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-ULP-1 Polyclonal Antibody
Target	ULP-1
Immunogen	Recombinant Yeast ULP 1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, ELISA, ICC/IF
Cellular Localization	Nuclear
Format	Liquid
Size	500 µg
Buffer	0.15M Sodium chloride, 0.02M Potassium phosphate. pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

ULP1, ubiquitin-like protein-specific protease 1, initially processes Smt3 and also acts as a deconjugating enzyme for Smt3 [Saccharomyces cerevisiae(Baker's yeast)]. Covalent modification of cellular proteins by the ubiquitin-like modifier SUMO (small ub

---

**Keywords**

Probable sentrin specific protease,Ubiquitin Like Protease,

---