



Rabbit Anti-Sheep Serum Polyclonal antibody (DPATB-H82812)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Sheep Serum Polyclonal Antibody
Target	Sheep Serum
Immunogen	pooled normal sheep serum
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Sheep
Purification	IgG fraction
Conjugate	Unconjugated
Applications	IEP
Format	Lyophilised: To one vial of lyophilized powder, add 3ml of deionized water. Rotate vial gently until powder dissolves.
Size	2 ml
Buffer	0.01M PBS, pH 7.2; preservative has been added.
Preservative	None
Storage	Store at +4°C.

BACKGROUND

Introduction serum collected from healthy normal sheep.

Keywords sheep serum,
